




MONAI ハンズオンセミナー

Colleen Ruan, Senior Solution Architect, Healthcare, Japan

Outline


- MONAIのご紹介（30分頃）
- MONAI Coreを利用して前処理、モデル学習を行う（1時間—2時間）
- MONAI Labelを利用してデータの自動アノテーションのデモンストレーション（30分頃）
- Q&A（15分）

MONAIのご紹介



NVIDIA CLARA


GENOMICS



NLP

Fever **PROBLEM** and urinary symptoms **PROBLEM** : A preliminary diagnosis of pyelonephritis **PROBLEM** was established. Other causes of fever **PROBLEM** were possible but less likely. The patient was **hypotensive PROBLEM** on initial assessment **TEST** with a blood pressure **TEST** of 80/40. Serum lactate **TEST** was **elevated PROBLEM** at 6.1. A bolus of IV fluid **TREATMENT** was administered (1.5L) but the patient remained **hypotensive PROBLEM**. Our colleagues from ICU were consulted. An arterial line **TREATMENT** was inserted for hemodynamic monitoring **TEST**. Hemodynamics were supported with levophed **TREATMENT** and crystalloids **TREATMENT**. Piptazo **TREATMENT** was started after blood and urine cultures **TEST** were drawn. After 12 hours serum lactate **TEST** had normalized and hemodynamics **TEST** had stabilized. Blood cultures **TEST** were positive for E Coli **PROBLEM** that was sensitive to all antibiotics **TREATMENT**. The patient was stepped down to oral ciprofloxacin **TREATMENT** to complete a total 14 day course of antibiotics **TREATMENT**. Fever **PROBLEM** and urinary symptoms **PROBLEM** : A preliminary diagnosis of pyelonephritis **PROBLEM** was established. Other causes of fever **PROBLEM** were possible but less likely. The patient was **hypotensive PROBLEM** on initial assessment **TEST** with a blood pressure **TEST** of 80/40. Serum lactate **TEST** was **elevated PROBLEM** at 6.1. A bolus of IV fluid **TREATMENT** was administered (1.5L) but the patient


IMAGING




INSTRUMENTS



CONVERSATIONAL AI




DRUG DISCOVERY




Medical Imaging is the essential instrument for Healthcare

2,000,000 Devices | 16,000 Companies | 10,000 Modalities




30% of World's Data is Healthcare | 90% of healthcare data is from medical imaging

MONAI: THE ESSENTIAL FRAMEWORK FOR THE MEDICAL AI ECOSYSTEM



MONAI Advisory Board

Bringing Together Domain Experts & Leaders



Stephen Aylward

Chair of the Advisory Board

Sebastien Ourselin

Klaus Maier-Hein

Jayashree Kalpathy-Cramer

Jorge Cardoso

Daniel Rubin

Lena Maier-Hein



Kevin Zhou

Nassir Navab

Andrew Feng

Nasir Rajpoot

Justin Kirby


Keyvan Farahani

Selnur Erdal

MONAI stands for Medical Open Network for AI and is a collaborative open-source projects for building the best practices for deep learning in Medical Imaging

MONAI

An AI Acceleration Framework for Medical Imaging



Frederick National Laboratory
for Cancer Research




NHS
Guy's and St Thomas'
NHS Foundation Trust



NIH
NATIONAL CANCER INSTITUTE



UNITED IMAGING



MONAI Impact

Standard for Imaging Research and Clinical Deployment

50K
Download per Month

450+
Github Projects

150+
Papers Published

2020

MONAI Core 0.1

Domain specific training library

4 publications

550 Bootcamp registrants

2022

MONAI 1.0

Model Zoo, Active Learning, Auto3D

Won 6 challenges

MONAI Tutorial, MICCIA (Sept 22)

2019

MONAI Inception
MICCAI Shenzhen

NVIDIA and KCL
initiate MONAI
community


2021

MONAI Label, MONAI Deploy
Expansion of imaging AI workflow


Won 5 challenges
709 Bootcamp registrants

MONAI 1.0


New Features Across MONAI




Radiology



Pathology




Endoscopy



MONAI MODEL ZOO NEW

A Hub for Pre-Trained Imaging AI Models



The screenshot shows the MONAI Model Zoo homepage. At the top, there are logos for Charité, King's College London, NVIDIA, Vanderbilt University, Warwick University, and Mayo Clinic. Below the header is a navigation bar with links: HOME, FRAMEWORKS, DOCS, RESOURCES, MODEL ZOO (which is highlighted in green), and GITHUB. The main content area features a large graphic of a neural network with green nodes and connections, surrounded by green circles with white plus signs. To the left of the graphic, the text reads: "MONAI Model Zoo hosts a collection of medical imaging models in the MONAI Bundle format. The MONAI Bundle format defines portable describes of deep learning models. A bundle includes the critical information necessary during a model development life cycle and allows users and programs to understand the purpose and usage of the models." The title "MONAI Model Zoo" is at the top left.

Features	Benefits
15 and growing pre-trained models across CT, MR, Pathology, Endoscopy.	Jumpstart training workflows
Packaged as MONAI Bundle for publication to Model Zoo for 1-click	Establish common standard for reproducible research & collaboration
	Broaden reach and impact of research
	Leverage the SOTA pre-trained models for downstream clinical tasks

MONAI BUNDLE NEW

Open standard for Model Definition & Package




```
Execute training:  
  
python -m monai.bundle run training \  
  --meta_file configs/metadata.json \  
  --config_file configs/train.json \  
  --logging_file configs/logging.conf  
  
Override the train config to execute multi-GPU training:  
  
torchrun --standalone --nnodes=1 --nproc_per_node=2 -m monai.bundle run training \  
  --meta_file configs/metadata.json \  
  --config_file "[configs/train.json,'configs/multi_gpu_train.json']" \  
  --logging_file configs/logging.conf  
  
Override the train config to execute evaluation with the trained model:  
  
python -m monai.bundle run evaluating \  
  --meta_file configs/metadata.json \  
  --config_file "[configs/train.json,'configs/evaluate.json']" \  
  --logging_file configs/logging.conf  
  
Execute inference:  
  
python -m monai.bundle run evaluating \  
  --meta_file configs/metadata.json \  
  --config_file configs/inference.json \  
  --logging_file configs/logging.conf  
  
Verify the metadata format:  
  
python -m monai.bundle verify_metadata --meta_file configs/metadata.json --filepath eval/schema.json  
  
Verify the data shape of network:  
  
python -m monai.bundle verify_net_in_out network_def --meta_file configs/metadata.json --config_file configs/inference.json  
  
Export checkpoint to TorchScript file:  
  
python -m monai.bundle ckpt_export network_def \  
  --filepath models/model.ts \  
  --ckpt_file models/model.pt \  
  --meta_file configs/metadata.json \  
  --config_file configs/inference.json
```

- MONAI Bundle is a self-contained model package with pre-trained weights and all meta data
- Build python workflows via structured configs
- Develop once and deploy anywhere
- Ease of use & flexibility to override & customize configs
- Hybrid programming with support for config to python conversion

MONAI LABEL

Intelligent Labeling & Learning Tool




- Now supports streaming modality (2D) Endoscopy **New**
- MONAI Label Server making application ready for deploy and serve
- Client Viewers: 3DSlicer, OHIF, PAIR, Digital Slide Archive, QuPath, CVAT **New**
- Data Store support DICOM & File archive
- REST APIs

インタラクティブなアノテーション機能 DEEPGROW 3D


General availability with ability to re-train the model

1. 3Dの次世代医用画像注釈モデル
 2. 高品質のセグメンテーションを提供するための前景と背景のクリック
 3. これまでに見られなかった画像のモデル予測が大幅に改善
 4. FoviaビューアSDKに統合され、商用の臨床統合を可能
 5. 3DSlicer、OHIFに統合され、オープンソースコミュニティを実現



MONAI LABEL: Active Learning

Reduce the Cost of Labeling Imaging Data



Features	Benefits
Support 3D Segmentation for Radiology & 2D Endoscopy <small>NEW</small>	2D Endoscopy benchmark showcase 75% reduction in need for labeled data
Charite University reduced labeling cost by 50% using Epistemic algorithm	Image selection algorithms (Aleatoric, Epistemic and FAS) intelligently chooses which data to learn from
Human-in-the loop to annotate the most difficult images for clinical inputs	Train models that can learn more from less data.
1-click Train from Viewer UI	Focus on annotations that provide the highest gain in model performance and address areas where the model has uncertainty
	Seamlessly initiate training with high value annotated data bringing clinicians into data science loop



Dr. Matthew Jolley
Pediatric Cardiologist and Anesthesiologist
Children's Hospital of Philadelphia

Case Study: Researcher at AMC

*“What MONAI provides is a fundamental, scalable and powerful architecture to accelerate that process such that **we can focus on our application in a specific (clinical) domain**, rather than the creation of that fundamental architecture.”*

(September 2022)

The Client

- Treats pediatric patients with congenital heart disease
- Aims to employ precision medicine from research to the OR
- Started building his own AI/ML to create patient-specific models for planning surgeries and interventions

The Problem

- No two patients and intervention are alike
- 3D Image segmentation to create patient-specific models was laborious and lengthy
- Building the homegrown infrastructure was mundane, distracting and taking many years

The Solution



- Interact with, visualize and extract all information in patient's heart valve images
- Rapidly create models that inform the design of the correct surgical intervention the first time
- Accelerated imaging AI workflow, from labeling data, to segmenting images, and achieving the model they want

Why MONAI

- Freed his team to focus on solving clinical applications, not the infra
- Provided the boiler plates to build, scale and deploy AI faster from cardiology other surgical domains.
- Open-source community to share ideas, labels, code and implementation
- Extensible MONAI framework that integrates into 3D slicer

MONAI CORE 1.0


Standard of Imaging Research Framework



- Robust, backward compatible API design **New**
- Supports new MetaTensor API **New** to store Torch Tensors and meta information.
- Modular domain specialized components: network arch; transforms; losses; metrics; data handlers etc.
- Learning workflows: Supervised; Self-supervised; Federated Learning; Network architecture search
- Up to 15x faster performance as compared to vanilla implementations

AUTO3D SEGMENTATION


Low-code Framework to Train High-quality 3D Segmentation Models



- Auto3D is a self-contained solution to build 3D Segmentation models with minimal user inputs
- Supports state of the art algorithms: SwinUNETR; DiNTS; SegResNet out of the box and ability to customize
- Built for high Accuracy & efficiency Auto3D maximizes the productivity of data science teams.
- Auto3D based solution with no human interaction amongst MICCAI challenges winning solutions: HEKCTOR22, Instance22, ISLES22

MONAI FEDERATED LEARNING NEW


Privacy Preserving Collaborative AI development



Features	Benefits
Federated Learning Client Algo APIs	Ability to define a MONAI client app that can run on any FL platform. Enabler for FL Toolkits interoperability.
NVIDIA FLARE Integration	Nvidia FLARE, the federated platform developed by Nvidia already integrated with MONAI FL Client APIs
MONAI Bundle compatibility	Allow for seamlessly extending any bundle from MONAI zoo to a federated paradigm

MONAI DEPLOY

End-to-End AI Clinical Workflow



Features	Benefits
MONAI App SDK	Build containerized application from pre-trained models in < 20 minutes
MONAI Informatics Gateway ^{NEW}	Connect hospital information systems like PACS, EHRs and viewers over DICOM, FHIR and HL7
MONAI Workflow Manager ^{NEW}	Maintain AI governance for Hospital IT



Jorge Cardoso
CTO at London AI Center for value-based healthcare

Case Study

*“What MONAI brings is the final layer, which is the way to scalably apply AI across the entire medical record **across every patient pathway.**”*

The Client

- Helps hospitals across all departments safely and effectively use and adopt AI
- Tasked with integrating trained models into hospital operations, IT systems and patient care

The Problem

- Hospital/patient information is stored in disparate systems and formats even for the same patient
- Integration of AI is difficult to do across patient pathways
- Lack of a domain-specific, easy-to-use tools, OS and software applications that helps developer productivity in the medical space

The Solution

- A platform that brings and scales data science into the clinic as insights for clinicians
- Allow clinicians to be part of the AI process, from labeling to how it is deployed and used for patient care
- Integrate all different kinds of data about a patient into a single AI platform

Why MONAI

- Open, standardized, inclusive, interoperable, scalable, improves accessibility – across every patient pathway worldwide
- Do not have to reinvent the wheel; symbiotic relationship between clinics and AI that improves patient-specific outcomes
- Foundational, pre-trained models in Model Zoo that can be used off the shelf and finetuned for specific applications

MONAI RESOURCES

Join the Growing Community Now !

- . MONAI Website: <https://monai.io/>
- . MONAI Slack: <https://forms.gle/QTxJq3hFictp31UM9>
- . MONAI Docs:
 - . MONAI Core: <https://docs.monai.io/en/stable/>
 - . MONAI Label: <https://docs.monai.io/projects/label/en/latest/index.html>
 - . MONAI Deploy App SDK: <https://docs.monai.io/projects/monai-deploy-app-sdk/en/latest/>
- . MONAI Github: <https://github.com/Project-MONAI>
 - . MONAI Core: <https://github.com/Project-MONAI/MONAI>
 - . MONAI Label: <https://github.com/Project-MONAI/MONAILabel>
 - . MONAI Deploy: <https://github.com/Project-MONAI/monai-deploy>
- . MONAI YouTube: <https://www.youtube.com/c/Project-MONAI>
 - . Overview Videos, Deep Dive Series, Bootcamp and other event recordings
- . MONAI Twitter: <https://twitter.com/ProjectMONAI>
 - . Follow for the latest announcements
- . MONAI Medium: <https://monai.medium.com/>
 - . Read about our latest releases and our upcoming research interview series



ACCESS MONAI EVERYWHERE


Nvidia Enterprise | Cloud | Launchpad



NVIDIA Certified & DGX
Data Science at enterprise scale



Integrated for scale and ease of use
across CSP



Experience MONAI on Nvidia Launchpad

The background of the slide features a dark, abstract design composed of numerous glowing green lines. These lines vary in thickness and intensity, creating a sense of depth and motion. They form a complex, organic-looking structure that resembles a neural network or a microscopic view of a biological tissue. The overall aesthetic is modern and scientific.

MONAI Coreを利用して前処理、モデル学習を行う

MONAI CORE

An AI Acceleration Framework for Medical Imaging

The screenshot shows the MONAI Core API Reference page. At the top, there's a navigation bar with links for "What's New", "Highlights", "API Reference" (which is highlighted in blue), "Installation Guide", "Development", and "More". Below the navigation bar is a search bar with placeholder text "Search the docs ...". To the right of the search bar is a keyboard shortcut "ctrl + K". The main content area has a title "API Reference" in large blue text. To the left of the main content is a sidebar titled "Section Navigation" with a list of categories: Applications, Auto3dseg, Federated Learning, Model Bundle, Transforms, Loss functions, Network architectures, Metrics, Optimizers, Data, Engines, Inference methods, Event handlers, Visualizations, and Utilities. The "Transforms" category is highlighted with a green border. At the bottom of the page, there are "Previous" and "Next" links. The "Previous" link points to "Research and Application", and the "Next" link points to "Applications".

What's New Highlights **API Reference** Installation Guide Development More ▾

Search the docs ... `ctrl + K`

Section Navigation

- Applications
- Auto3dseg
- Federated Learning
- Model Bundle
- Transforms
- Loss functions
- Network architectures
- Metrics
- Optimizers
- Data
- Engines
- Inference methods
- Event handlers
- Visualizations
- Utilities

API Reference

Applications
Auto3dseg
Federated Learning
Model Bundle
Transforms
Loss functions
Network architectures
Metrics
Optimizers
Data
Engines
Inference methods
Event handlers
Visualizations
Utilities

Previous [Research and Application](#) Next [Applications](#)

MONAI COREハンズオンを開始!

https://github.com/colleenrpy/MONAI_hands-on

The screenshot shows a GitHub repository page for 'MONAI_hands-on'. At the top, there are navigation buttons for 'main' (selected), '1 branch', '0 tags', 'Go to file', 'Add file', and 'Code'. Below this is a list of files and folders:

File/Folder	Description	Last Commit
colleenrpy	Created using Colaboratory	5cc4dcf 2 hours ago
3d_segmentation	Created using Colaboratory	2 hours ago
MONAICore_v0.8.1	Add files via upload	4 days ago
MONAILabel_v0.4.dev2215/apps/pat...	Add files via upload	3 days ago
data	Add files via upload	4 days ago
End-To-End Workflow with MONAI.ip...	Created using Colaboratory	2 hours ago
Env.ipynb	Add files via upload	4 days ago
MONAI Datasets and Caching_soluti...	Created using Colaboratory	2 hours ago
MONAI_Transforms_Tutorial.ipynb	Created using Colaboratory	2 hours ago

At the bottom, there is a light blue box with the text 'Help people interested in this repository understand your project by adding a README.' and a green 'Add a README' button.

On the left side of the main content area, there is a text overlay 'まずはここをクリック' with a green arrow pointing towards the 'Add a README' button.

<https://docs.monai.io/en/stable/api.html>







MONAI Labelを利用してデータの自動アノテーションの
デモンストレーション

What is MONAI Label?


- An intelligent **open source** image labeling and learning tool that enables users to create annotated datasets and build AI annotation models for clinical evaluation
- Framework for developing and deploying **MONAI Label Apps to train and do inference using Deep Learning models**
- MONAI Label is one of the first frameworks that introduces **Active learning strategies in a software annotation setting**
- It is all Python and can be installed with simple “**pip install monailabel**”
- Supported viewers:
 - **3D Slicer**
 - **Open Health Imaging Foundation (OHIF)**



MONAI Label Infrastructure: server-client system



MONAI Label App: *Active Learning Strategies*




- **Active learning** is a semi-supervised machine learning approach where the algorithm can choose which data it wants to learn from
- Initially, random images are selected to be labeled

However, random is not always the most efficient.
- After having a pretrained model, uncertainty of each image can be computed. Unlabeled samples that need more attention from the clinician will be selected
- Selection of harder samples or samples that need more attention

Why use MONAI Label?

- **Researcher Perspective:** MONAI Label allows researchers to
 - Create new annotation methods
 - Rapid App prototyping
 - Implement active learning techniques
 - Verify their effectiveness in real-world scenarios
 - Make incremental improvements
 - Readily deploy labeling apps to wider audiences
- **Clinician Perspective:** MONAI Label reduces the time and effort of annotating new datasets
 - Ready-to-use **3DSlicer**
 - Pre-built plugin for **OHIF Viewer**

How to create a MONAI Label App?



MONAI Label Env Setup

Server

Hardware

- GPU: Pascal or newer, RAM >= 16Gb

Software

- PyTorch Docker Container (NGC):
<https://catalog.ngc.nvidia.com/orgs/nvidia/containers/pytorch>

Steps:

- Install MONAI Label by pip (MONAI will be installed automatically)
`$ pip install monailabel`
- Create your own app or download sample app by:
`$ monailabel apps --name radiology --download --output .`
- Create a dataset folder with your unlabeled data.
- Start the MONAI Label server with sample app:
`$ monailabel start_server --app radiology --studies /path/to/your/dataset --conf models deepeedit`

Client (3D Slicer)

<https://download.slicer.org/>

- 3D Slicer Installation



Windows



macOS




Linux

prerequisites

Stable Release	version 4.11.20210226 revision 29738 built 2021-03-01	version 4.11.20210226 revision 29738 built 2021-02-28	version 4.11.20210226 revision 29738 built 2021-02-27
Preview Release	version 4.13.0 revision 30785 built 2022-04-18	version 4.13.0 revision 30785 built 2022-04-18	version 4.13.0 revision 30785 built 2022-04-18

- MONAI Label Plugin Installation



- Connect to the MONAI Label server:
Enter the server URL including “http://”
<http://your.server.ip.address:port>

MONAI Label server:




MONAI Label Pre-trained Models

The default pre-trained model in the sample app is for **spleen segmentation**. If you want to try more pre-trained models, please follow the steps below:


1. Select a pre-trained segmentation model from NGC:

<https://catalog.ngc.nvidia.com/models?filters=&orderBy=dateModifiedDESC&query=clara>



2. Copy the pre-trained weights into MONAI Label APP's model folder.
3. According to the configs of the pre-trained model, please update the APP's trainer, infer and config.

MONAI Label With Pre-trained Models



MONAI Label Without Pre-trained Models

If we don't have a pre-trained model before annotation, we can also apply MONAI Label to accelerate the annotating process.

1. Create your **custom MONAI Label app**. You should define the **model architecture** you want to use and the **appropriate transforms**. (Including pre-processing, augmentation and post-processing), then make the model to be initialized randomly.
2. Annotate a few data by yourself and submit to MONAI Label server, then train the initialized model by the annotated data.
3. Now you have a pre-trained model. Continuing looping the Active Learning Strategy to annotate your dataset.

